

## Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {HRP} (STJS000129)

STJS000129

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	HRP-conjugated rabbit polyclonal anti-Guinea Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA/WB/IHC/ICC
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Guinea Pig

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	HRP
<b>Purification</b>	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
<b>Formulation</b>	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>