

Mouse Anti-Chicken IgY Heavy & Light Chain antibody {HRP} (STJS000140)

STJS000140

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated mouse monoclonal anti-Chicken IgY Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA/WB/IHC/ICC
Host/Source	Mouse
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	HRP
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgY
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence