

Rabbit Anti-Sheep IgG Heavy & Light Chain antibody {Biotin} (STJS000171)

STJS000171

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated rabbit polyclonal anti-Sheep IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA/WB/IHC/ICC
Host/Source	Rabbit
Reactivity	Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Biotin
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence