

**Rabbit Anti-Human IgM antibody {FITC} (STJS000198)**

STJS000198

**GENERAL INFORMATION**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	FITC-conjugated rabbit polyclonal anti-Human IgM secondary antibody. For use in most research applications.
<b>Applications</b>	IF/FC
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	FITC
<b>Purification</b>	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IF 1:250-1:2000 FCM 1:250-1:2000
<b>Formulation</b>	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen Region**  
**Specificity**  
**Immunogen Sequence**