

Rabbit Anti-Monkey IgG Heavy & Light Chain antibody {FITC} (STJS000199)

STJS000199

GENERAL INFORMATION

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Secondary antibodies |
| Short Description | FITC-conjugated rabbit polyclonal anti-Monkey IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| Applications | IF/FC |
| Host/Source | Rabbit |
| Reactivity | Monkey |

PRODUCT PROPERTIES

| | |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | FITC |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range | IF 1:250-1:2000 FCM 1:250-1:2000 |
| Formulation | Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |