

Goat Anti-Mouse IgG Highly Cross-Adsorbed (H&L) antibody {Alexa Fluor 633} (STJS000838)

STJS000838

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 633-conjugated goat polyclonal anti-Mouse IgG Highly Cross-Adsorbed H&L secondary antibody. For use in most research applications.
Applications	
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Alexa Fluor 633
Purification	Affinity purified
Dilution Range	1-10 Mu g/mL of the IgG conjugate for most applications. Appropriate dilutions of the conjugate should be determined empirically.
Formulation	Blue solution in pH~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.
Isotype	IgG
Storage	Store at 4°C for several months. Protect from light. For longer storage, divide the conjugate into small aliquots and freeze at -20°C.
Instruction	Avoid repeated freezing and thawing.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	H&L
Immunogen Region	H&L
Specificity	
Immunogen Sequence	