

Rabbit Anti-Guinea Pig IgG H&L antibody {AMCA} (STJS000881)

STJS000881

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	AMCA-conjugated rabbit polyclonal anti-Guinea Pig IgG H&L secondary antibody. For use in most research applications.
Applications	ELISA/IHC/IF/FC
Host/Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	AMCA
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Guinea Pig IgG
Immunogen Region	
Specificity	
Immunogen Sequence	