

Goat Anti-Mouse IgM mu chain antibody {TRITC} (STJS000914)

STJS000914

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	TRITC-conjugated goat polyclonal anti-Mouse IgM mu chain secondary antibody. For use in most research applications.
Applications	ELISA/IHC/IF/FC
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	TRITC
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Mouse IgM mu chain
Immunogen Region	
Specificity	
Immunogen Sequence	