

**Rabbit Anti-Human IgG Gamma chain antibody {TRITC} (STJS000928)**

STJS000928

**GENERAL INFORMATION**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	TRITC-conjugated rabbit polyclonal anti-Human IgG Gamma chain secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA/IHC/IF/FC
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	TRITC
<b>Purification</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
<b>Formulation</b>	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Human IgG gamma chain
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	