

Rabbit Anti-Chicken IgY H&L antibody {Cy3} (STJS001005)

STJS001005

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Cy3-conjugated rabbit polyclonal anti-Chicken IgY H&L secondary antibody. For use in most research applications.
Applications	ELISA/IHC/IF/FC
Host/Source	Rabbit
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Cy3
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Chicken IgY
Immunogen Region	
Specificity	
Immunogen Sequence	