

**Rabbit Anti-Cow IgG H&L antibody {Cy5} (STJS001046)**

STJS001046

**GENERAL INFORMATION**

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Secondary antibodies   |
| <b>Short Description</b> | Cy5-conjugated rabbit polyclonal anti-Cow IgG H&L secondary antibody. For use in most research applications. |
| <b>Applications</b>      | ELISA/IHC/IF/FC  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Cow  |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         | Cy5   |
| <b>Purification</b>        | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| <b>Dilution Range</b>      | IHC 1:200-1:1000<br>IF 1:200-1:1000<br>FCM 1:2000-1:4000<br>ELISA   |
| <b>Formulation</b>         | Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide  |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                    |

**TARGET INFORMATION**

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            |            |
| <b>Gene Symbol</b>        |            |
| <b>Uniprot ID</b>         |            |
| <b>Immunogen</b>          | Bovine IgG |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |