

Rabbit Anti-Human IgG Gamma chain antibody {Cy5.5} (STJS001100)

STJS001100

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Cy5.5-conjugated rabbit polyclonal anti-Human IgG Gamma chain secondary antibody. For use in most research applications.
Applications	ELISA/IHC/IF/FC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Cy5.5
Purification	The antibody was isolated from antiserum by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000
	IF 1:200-1:1000
	FCM 1:2000-1:4000
	ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human IgG gamma chain
Immunogen Region	
Specificity	
Immunogen	
Sequence	