

Mouse Anti-Monkey IgG H&L antibody {Cy5.5} (STJS001115)

STJS001115

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Cy5.5-conjugated mouse monoclonal anti-Monkey IgG H&L secondary antibody. For use in most research applications.
Applications	ELISA/IHC/IF/FC
Host/Source	Mouse
Reactivity	Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	Cy5.5
Conjugation	
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Monkey IgG
Immunogen Region	
Specificity	
Immunogen Sequence	