

Goat Anti-Mouse IgG H+L antibody {AP} (STJS001171)

STJS001171

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	
Applications	WB
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	AP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:50000
Formulation	PBS with 0.025% Sodium Azide, 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Mouse IgG
Immunogen Sequence	