

Rabbit Anti-Goat IgG H+L antibody {HRP} (STJS001178)

STJS001178

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	
Applications	WB/IHC-P
Host/Source	Rabbit
Reactivity	Goat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:100000 IHC-P:1:500-1:5000
Formulation	PBS with 0.025% Sodium Azide, 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Goat IgG
Immunogen Sequence	