

Goat Anti-IgG2a antibody {HRP} (STJS001201)
STJS001201

GENERAL INFORMATION


Product Type	Secondary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Goat
Reactivity	

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:2000-1:5000 ELISA:1:2000-1:5000
Formulation	PBS with 50% Glycerol, pH 7.4.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Mouse IgG2a
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081