

**Rabbit Anti-IgM Fc antibody {HRP} (STJS001217)**  
STJS001217

**GENERAL INFORMATION**


<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SRM
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	HRP
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:10000-1:50000
<b>Formulation</b>	PBS with 50%Glycerol, 0.2%BSA, 0.03%Proclin300, Negative IgG, pH 7.3.
<b>Isotype</b>	
<b>Storage Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	Mouse IgM
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081