

Anti-ASCC3 antibody (1-111) (STJ110267)

STJ110267

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASCC3 (1-111) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10973
Gene Symbol	ASCC3
Uniprot ID	ASCC3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ASCC3 (NP_071374.1).
Immunogen Region	1-111
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of A549 cells using ASCC3 antibody (STJ110267). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081