

Anti-RNF125 antibody (1-232) (STJ110289)

STJ110289

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RNF125 (1-232) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54941
Gene Symbol	RNF125
Uniprot ID	RN125_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human RNF125 (NP_060301.2).
Immunogen Region	1-232
Specificity	
Immunogen Sequence	