

Anti-LAMB3 antibody (18-230) (STJ110426)

STJ110426

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LAMB3 (18-230) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3914
Gene Symbol	LAMB3
Uniprot ID	LAMB3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-230 of human LAMB3 (NP_000219.2).
Immunogen Region	18-230
Specificity	
Immunogen Sequence	