

Anti-WNT5B antibody (120-359) (STJ110611)

STJ110611

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WNT5B (120-359) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	81029
Gene Symbol	WNT5B
Uniprot ID	WNT5B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-359 of human WNT5B (NP_116031.1).
Immunogen Region	120-359
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using WNT5B antibody (STJ110611) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.



Immunohistochemistry of paraffin-embedded rat lung using WNT5B antibody (STJ110611) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat bone marrow using WNT5B antibody (STJ110611) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat heart using WNT5B antibody (STJ110611) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081