

Anti-KCNJ4 antibody (316-445) (STJ110732)

STJ110732

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KCNJ4 (316-445) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3761
Gene Symbol	KCNJ4
Uniprot ID	KCNJ4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 316-445 of human KCNJ4 (NP_690607.1).
Immunogen Region	316-445
Specificity	
Immunogen Sequence	