

Anti-STX5 antibody (55-355) (STJ110753)

STJ110753

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-STX5 (55-355) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6811
Gene Symbol	STX5
Uniprot ID	STX5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-355 of human STX5 (NP_003155.2).
Immunogen Region	55-355
Specificity	
Immunogen Sequence	