

Anti-TFAP2A antibody (1-437) (STJ110952)

STJ110952

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TFAP2A (1-437) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:1000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage	Store in a freezer at -20°C and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7020
Gene Symbol	TFAP2A
Uniprot ID	AP2A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-437 of human Activator protein 2 (AP-2/Activator protein 2 (AP-2/TFAP2A)) (NP_001027451.1).
Immunogen Region	1-437
Specificity	
Immunogen Sequence	