

Anti-OBSCN antibody (STJ11100081)

STJ11100081

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OBSCN is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	84033
Gene Symbol	OBSCN
Uniprot ID	OBSCN_HUMAN
Immunogen	Recombinant protein of human OBSCN.
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U-2 OS cells using OBSCN antibody (STJ11100081) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of rat heart using OBSCN antibody (STJ11100081) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse heart using OBSCN antibody (STJ11100081) at dilution of 1:100. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081