

Anti-NUDCD3 antibody (STJ11100390)

STJ11100390

GENERAL INFORMATION

| | |
|--------------------------|-----------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-NUDCD3 is suitable for use in Western Blot. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|-----------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------------|
| Gene ID | 23386 |
| Gene Symbol | NUDCD3 |
| Uniprot ID | NUDC3_HUMAN |
| Immunogen | Recombinant protein of human NUDCD3. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |