

Anti-PLEKHH3 antibody (STJ11100452)

STJ11100452

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PLEKHH3 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79990
Gene Symbol	PLEKHH3
Uniprot ID	PKHH3_HUMAN
Immunogen	Recombinant protein of human PLEKHH3.
Immunogen Region	
Specificity	
Immunogen Sequence	