

Anti-OSBPL9 antibody (STJ11100480)

STJ11100480

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OSBPL9 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114883
Gene Symbol	OSBPL9
Uniprot ID	OSBL9_HUMAN
Immunogen	Recombinant protein of human OSBPL9.
Immunogen Region	
Specificity	
Immunogen Sequence	