

Anti-C14orf37 antibody (700-774) (STJ11100498)

STJ11100498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C14orf37 (700-774) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	145407
Gene Symbol	ARMH4
Uniprot ID	ARMD4_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-774 of human C14orf37 (NP_001001872.2).
Immunogen Region	700-774
Specificity	
Immunogen Sequence	