

Anti-CXCL5 antibody (37-114) (STJ11100886)

STJ11100886

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CXCL5 (37-114) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6374
Gene Symbol	CXCL5
Uniprot ID	CXCL5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-114 of human CXCL5 (NP_002985.1).
Immunogen Region	37-114
Specificity	
Immunogen Sequence	