

Anti-MAS1 antibody (246-325) (STJ11100947)

STJ11100947

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAS1 (246-325) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4142
Gene Symbol	MAS1
Uniprot ID	MAS_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 246-325 of human MAS1 (NP_002368.1).
Immunogen Region	246-325
Specificity	
Immunogen Sequence	