

Anti-Phospho-CTNND1-S320 antibody (STJ11101117)

STJ11101117

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Delta-Catenin/p120 Catenin-S320 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1500
Gene Symbol	CTNND1
Uniprot ID	CTND1_HUMAN
Immunogen	A phospho synthetic peptide corresponding to residues surrounding S320 of human Delta-Catenin/p120 Catenin.
Immunogen Region	
Specificity	
Immunogen Sequence	