

## Anti-Phospho-STK3/STK4-T180/T183 antibody (STJ11101128)

STJ11101128

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-STK3/STK4-T180/T183 is suitable for use in Western Blot.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000
<b>Formulation</b>	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store in a freezer at -20°C and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">6789</a> <a href="#">6788</a>
<b>Gene Symbol</b>	<a href="#">STK4</a> <a href="#">STK3</a>
<b>Uniprot ID</b>	<a href="#">STK4_HUMAN</a> <a href="#">STK3_HUMAN</a>
<b>Immunogen</b>	A phospho specific peptide corresponding to residues surrounding T180/T183 of human STK3/STK4.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot analysis of extracts of various cell lines, using Phospho-STK3/STK4-T180/T183 rabbit polyclonal antibody (STJ11101128) at 1:1000 dilution. Both NIH/3T3 cells and C6 cells were treated by Staurosporine (1 $\mu$ M) at room temperature for 3 hours. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 $\mu$ g per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081