

Anti-LOC101736823 antibody (251-460) (STJ11101474)

STJ11101474

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LOC101736823 (251-460) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	B.mori

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 251-460 of Bombyx mori LOC101736823.
Immunogen Region	251-460
Specificity	
Immunogen Sequence	