

Anti-MGST3 antibody (1-100) (STJ11101502)

STJ11101502

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-MGST3 (1-100) is suitable for use in Western Blot. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 4259 |
| Gene Symbol | MGST3 |
| Uniprot ID | MGST3_HUMAN |
| Immunogen | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MGST3 (NP_004519.1). |
| Immunogen Region | 1-100 |
| Specificity | |
| Immunogen Sequence | |