

Anti-MAPK8/MAPK9/MAPK10 antibody [ARC0291] (STJ11102174)

STJ11102174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-JNK1/2/3 is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	ARC0291
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID [5599](#)
[5602](#)
[5601](#)
[MAPK8](#)
[MAPK10](#)
[MK08_HUMAN](#)
[MK10_HUMAN](#)
[MK09_HUMAN](#)
[<](#)

Immunogen A synthesized peptide derived from human JNK1/2/3

Immunogen Region
Specificity
Immunogen Sequence



Western blot analysis of extracts of various cell lines, using JNK1/2/3 rabbit monoclonal antibody (STJ11102174) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence analysis of mouse brain using JNK1/2/3 rabbit monoclonal antibody (STJ11102174) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081