

Anti-Phospho-Phospho-PIK3C3-S165-S165 antibody (STJ11102669)

STJ11102669

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PIK3C3-S165 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5289
Gene Symbol	PIK3C3
Uniprot ID	PK3C3_HUMAN
Immunogen	A synthetic phosphorylated peptide around S165 of human PIK3C3/VPS34 (NP_002638.2).
Immunogen Region	
Specificity	
Immunogen Sequence	