

Anti-PPP1R15A antibody (STJ11102850)

STJ11102850

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PPP1R15A is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23645
Gene Symbol	PPP1R15A
Uniprot ID	PR15A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 375-674 of human PPP1R15A (NP_055145.3).
Immunogen Region	
Specificity	
Immunogen Sequence	