

Anti-TAS2R14 antibody (STJ11102859)

STJ11102859

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAS2R14 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50840
Gene Symbol	TAS2R14
Uniprot ID	T2R14_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human TAS2R14 (NP_076411.1).
Immunogen Region	
Specificity	
Immunogen Sequence	