

Anti-WWC3 antibody (STJ11102905)

STJ11102905

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WWC3 is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, pH7.3.
Isotype	IgG
Storage Instruction	Store at 2-8°C

TARGET INFORMATION

Gene ID	55841
Gene Symbol	WWC3
Uniprot ID	WWC3_HUMAN
Immunogen	A synthetic peptide of human WWC3.
Immunogen Region	
Specificity	
Immunogen Sequence	