

Anti-NPY5R antibody (STJ11103226)

STJ11103226

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPY5R is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4889
Gene Symbol	NPY5R
Uniprot ID	NPY5R_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NPY5R (NP_001304020.1).
Immunogen Region	
Specificity	
Immunogen Sequence	