

Anti-ANKRD27 antibody (STJ11103248)

STJ11103248

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANKRD27 is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84079
Gene Symbol	ANKRD27
Uniprot ID	ANR27_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 560-710 of human ANKRD27 (NP_115515.2).
Immunogen Region	
Specificity	
Immunogen Sequence	