

Anti-AMER1 antibody (STJ11103307)

STJ11103307

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-AMER1 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	139285
Gene Symbol	AMER1
Uniprot ID	AMER1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-660 of human AMER1 (NP_689637.3).
Immunogen Region	
Specificity	
Immunogen Sequence	