

Anti-PIK3IP1 antibody (STJ11103409)

STJ11103409

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PIK3IP1 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:1000
Formulation	PBS containing 0.02% Sodium Azide, pH7.3.
Isotype	IgG
Storage Instruction	Store at 2-8°C

TARGET INFORMATION

Gene ID	113791
Gene Symbol	PIK3IP1
Uniprot ID	P3IP1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human PIK3IP1 (NP_443112.2).
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of SH-SY5Y cells, using PIK3IP1 antibody (STJ11103409) at 1:500 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081