

Anti-PDS5A antibody (STJ11103567)

STJ11103567

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDS5A is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:1000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23244
Gene Symbol	PDS5A
Uniprot ID	PDS5A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1120 of human PDS5A (Q29RF7).
Immunogen Region	
Specificity	
Immunogen Sequence	