

Anti-MRPL37 antibody (STJ11103575)

STJ11103575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MRPL37 is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51253
Gene Symbol	MRPL37
Uniprot ID	RM37_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 116-423 of human MRPL37 (NP_057575.2).
Immunogen Region	
Specificity	
Immunogen Sequence	