

Anti-S100/S100A1 antibody (STJ11103624)

STJ11103624

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-S100/S100A1 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6271
Gene Symbol	S100A1
Uniprot ID	S10A1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human S100/S100/S100A1 (NP_006262.1).
Immunogen Region	
Specificity	
Immunogen Sequence	