

Anti-GST-Tag antibody {HRP} (STJ11103703)

STJ11103703

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-GST-Tag is suitable for use in Western Blot.
Applications	WB
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB 1:2000-1:10000
Formulation	PBS containing 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-218 of GST protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot analysis of recombinant GST protein using HRP-conjugated mouse anti-GST-Tag antibody (STJ11103703) at dilution of 1:5000. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081