

Anti-Phospho-VANGL2-S11/S14 antibody (STJ11103740)

STJ11103740

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Phospho-VANGL2-S11/S14 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57216
Gene Symbol	VANGL2
Uniprot ID	VANGL2_HUMAN
Immunogen	A phospho specific peptide corresponding to residues surrounding S11/S14 of human VANGL2
Immunogen Region	
Specificity	
Immunogen Sequence	