

## Anti-Phospho-AKT1-S308 & S309 & S305 antibody (STJ11103779)

STJ11103779

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-AKT1-S308 & S309 & S305 is suitable for use in Western Blot.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000
<b>Formulation</b>	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store in a freezer at -20°C and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">208</a>
	<a href="#">207</a>
	<a href="#">10000</a>
	<a href="#">AKT2</a>
	<a href="#">AKT1</a>
	<a href="#">AKT2_HUMAN</a>
	<a href="#">AKT1_HUMAN</a>
	<a href="#">AKT3_HUMAN</a>
	<

**Immunogen**  
**Immunogen Region**  
**Specificity**  
**Immunogen Sequence**

A phospho synthetic peptide corresponding to residues surrounding S308/S309/S305 of human AKT1/AKT2/AKT3.



Western blot analysis of extracts of Jurkat cells, using Phospho-Akt-T308 antibody (STJ11103779) at 1:1000 dilution. Jurkat cells were treated by Serum-starvation overnight at 37 °C. Jurkat cells were treated by Calyculin A (100 nM) at 37 °C for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081